

etac[®]
Creating Possibilities



etac[®]

Etac Socky

Stocking aid

Socky

Getting dressed should be simple, but pain or limited mobility can make putting on socks or stockings difficult. The Socky stocking aids are specifically designed to assist with this challenge, offering three easy-to-use models, each with slightly different functions. These practical aids provide an effective solution for maintaining independence while dressing.

Socky short

Socky short is a practical stocking aid for both short and knee long socks including flight socks. It consists of a plastic base with two hooks under a pre-folded nylon fabric, making it easier to slide socks on.



Socky long

Ideal for those with limited hand strength or mobility, Socky long reduces friction between their socks and the floor or carpet. It consists of a long plastic sheet beneath a pre-folded nylon fabric, and includes a convenient handle on top. Socky long is suitable for both short and knee-long socks.




Socky support


This stocking aid is designed to suit different lengths and types of support and compression stockings. Socky support consists of two long nylon fabric sections that connect with hook-and-loop fasteners. Each part has a large handle to grip to pull the stocking over the foot and leg.



Socky	Item no.	Weight
Socky short	80601002-2	Socky short: 94 gr
Socky long	80601001-2	Socky long: 146 gr
Socky support	80601003-2	Socky support: 90 gr

Length	Maintenance
Socky short: 400 mm	 Handwash the product in max. temperature 40°C
Socky long: 600 mm	
Socky support: 1000 mm	

Materials
Fabric: Polyester
Plastic parts: Polypropylene



Follow these steps for a perfect result



1. Pull up the sock properly and over the hooks



2. Place the string behind the calf



3. Pull up with ease





Etac is a world-leading developer of innovative assistive devices and patient handling equipment. Since 1973, we have been committed to improving quality of life for the individual, family members and caregivers.

For the intended purpose of the products and latest information, please visit www.etac.com

Etac

- +46 371 587 330
- info@etac.se
- www.etac.com

